

Analysis of Isothermal Crystallization Process of High Polymers with DSC

Crystalline high polymers including polyethylene, polypropylene and nylon produce crystallized structures if they are left alone for a long time at a given constant temperature lower than that when they have been completely fused. This phenomenon is called isothermal crystallization.

Measurement of the isothermal crystallization process serves to determine a half-life ($t_{1/2}$), which can be used as a yardstick for crystallization speed, evaluation of core substance of crystallization, and the study of bridging effect.

DSC, excellent in thermal sensibility, is an indispensable piece of equipment for measuring the process of isothermal crystallization where a specimen is rapidly cooled after fusion and kept under a strictly controlled temperature.

The following is a report on measurements of isothermal crystallization process of polyethylene and polypropylene with DSC.

A typical process of isothermal crystallization is presented in Fig.1. The curve shows that a fused specimen is rapidly cooled down to a given temperature, kept at a constant temperature, and produces an exothermic peak of crystallization after a lag time "tid."

ΔH_c of whole crystallization calories can be known from the whole peak area, and crystallization calories at a given time t are calculated by the following equation.

$$1-\theta = \frac{\int_{t_{id}}^t \frac{dH}{dt} \cdot dt}{\Delta H_c}$$

where θ : Ratio of non-crystallized portion,

ΔH_c : whole crystallization calories.

Values of q for each point of t make an isothermal crystallization curve that can be drawn from the values θ determined by above equation for each t .







Fig.1
A Typical Curve of Isothermal Crystallization

Analysis of isothermal crystallization process of polypropylene

Figs.2-4 show results of analysis of isothermal crystallization process of polypropylene held at different temperatures. Each specimen weighing 10.62mg was rapidly cooled at a rate of $-50^{\circ}\text{C}/\text{min}$ after fusion and held at constant temperature.






Fig. 5
Isothermal Crystallization Curves of Polypropylene

Analysis of polyethylene

Similarly to polypropylene, each specimen weighing 10.84mg was rapidly cooled at a rate of $-50^{\circ}\text{C}/\text{min}$ after fusion, and held at a constant temperature.




Fig.6 Curve for 122°C




Fig.7 Curve for 124°C

Fig.8 shows isothermal crystallization curves of polyethylene.




Fig.8 Isothermal Curves of Polyethylene

* Please be advised that data obtained before the implementation of the current Weights and Measures Law may be presented in terms of gravimetric unit.



SHIMADZU CORPORATION International Marketing Division
3,Kanda-Nishikicho 1-chome, Chiyoda-ku,Tokyo 101-8448,
JapanPhone: 81(3)3219-5641 Fax: 81(3)3219-5710
Cable Add,;SHIMADZU TOKYO

SHIMADZU SCIENTIFIC INSTRUMENTS, INC.
7102, Riverwood Drive, Columbia, Maryland 21046, U.S.A.
Phone: 1(410)381-1227 Fax: 1(410)381-1222 Toll Free:1(800)477-1227

SHIMADZU DEUTSCHLAND GmbH

Albert-Hahn-Strasse 6-10, D-47269 Duisburg, F.R. Germany
Phone: 49(203)7687-0 Fax: 49(203)7666-25

SHIMADZU (ASIA PACIFIC) PTE LTD.
16 Science Park Drive #01-01, Singapore Science Park, Singapore 118227, Republic of Singapore
Phone: 65-778-6280 Fax: 65-779-2935

SHIMADZU SCIENTIFIC INSTRUMENTS(OCEANIA)PTY. LTD.
Rydalmere Business Park, Unit T, 10-16, South Street Rydalmere, N.S.W. 2116, Australia
Phone: 61(2)9684-4200 Fax: 61(2)9684-4055

SHIMADZU DO BRASIL COMERCIO LTDA.
Rua Cenzo Sbrighi, 25 CEP 05036-010 Agua Branca, Sao Paulo, BRAZIL
Phone: 55(11)3611-1688 Fax: 55(11)3611-2209

SHIMADZU (HONG KONG) LIMITED.
Suite 1028, Ocean Centre, Harbour City, Tsim Sha tsui, Kowloon, HONG KONG
Phone: 852(2375)4979 Fax: 852(2199)7438

Overseas Office
Istanbul, Beijing, Shanghai, Guangzhou, Shenyang, Chengdu, Moscow